

Solid Xenon

Jonghee Yoo



**Fermilab Center for Particle Astrophysics Retreat
18 April 2009**

CDMS Axion Search

2




See arXiv:0902.4693 (submitted to PRL)





- Searching for the invisible axion models
- For discovery or solar model consistent search
 $\Rightarrow g_{a\gamma\gamma} < 0.5 \times 10^{-10} \text{ GeV}^{-1}$ (solar ν limit)
- Require 100 kg of detector fiducial volume
 and $< 10^{-4}$ dru of gamma/e background level
 (dru = counts/kg/day/keV)
- The current CDMS gamma background level: ~ 1 dru

Big Picture : Particle & Astro

3





- Most models are at $\sigma > 10^{-47} \text{cm}^2$
- The goal for the major Dark Matter experiments:
 - **nuclear recoil background** level should be controlled
less than $\sim 10^{-7} \text{ dru}$ (= counts/day/kg/keV)
 - Current CDMS **nuclear recoil background level** (world best) : $\sim 10^{-4} \text{ dru}$
- 10^4 times better background control required (with a huge detector volume)!

10^{-7} counts/day/kg/keV

1 counts/3,000years/kg/10keVs

1 counts/year/3tons/10keVs


- When you first build a detector (with raw target material and signal readout)
you will typically have 10~100 Hz of low energy event rate
by cosmic rays, ambient gammas and internal radioactive contaminations etc
- In order to realize a **low background detector** (1 event/years/tons/keVs ultra-silent detector)
 - $\sim 10^8$ background reduction from the raw detector
 - Gas distillation/purification, material screening,
deeper site, active/passive shielding ...
 - **Extremely Challenging !**
- **Self-shielding** is suggested as a solution
 - Shield out gammas and neutrons at the surface of the detector target material
 - Noble liquid detectors : XENON / LUX / XMASS / WArP / ArDM / MAX / ...
- Concerns about the noble liquid detectors
 - ☑ Background contamination from PMTs ? → QUPID, acrylic vessels ...
 - ☑ Electron drift above meter scale? → electro-negative purity control ...
 - ☑ 2nd contamination from the circulation loop ? → closed circulation system ...
 - ☑ Outgas or vacuum leak? → do better job

Fine Tuning the Existing Recipes

7



100 kg Dual-phase Xenon TPC
Waseda University (~2002?)
Kamioka@Japan



170 kg Dual-phase Xenon TPC
XENON Collaboration (~2009)
Gran Sasso@Italy




300 kg Dual-phase Xenon TPC
LUX Collaboration (~2009)
SUSEL@USA

Fine Tuning the Existing Recipes

8

XMASS : **Single phase** detector (scintillation light readout)



- Kamioka Underground Laboratory, Japan
- 100 kg R&D phase was successful
 - Position and energy, self shielding, BG study
- Background reduction plan
 - PMT gammas : self shielding, fiducial volume cut
10 mBq/PMT (Hamamatsu R10789) achieved!
 - External gammas and neutrons : reduce with water shielding
- 800 kg (100 kg FV) operation will start in 2009



100 kg R&D detector



800 kg DM detector




Two Fundamental Questions

9

1. How will you **make** such a low background detector ?
2. How will you **keep** the detector pure ?

Solid Xenon



Why Xenon ?


- No long-lived Xe radio isotope (no intrinsic background)
- High yield of scintillation light
- Scintillation wavelength : 175nm (optically transparent)
- Relatively high melting point : $T_m = 161\text{K}$
- Simple crystal structure : fcc (same with Ge)
- Easy purification (distillation, etc)
- Self shielding : $Z=54$

Why Solid ?

- For solar axion search, being a crystal is crucial (Bragg scattering)
- Even more scintillation light (61 γ / keV) than LXe (42 γ / keV)
- Drifting electrons is easier in the crystal
- Superb superconducting sensors are running at low temperature (mK ~ K)
- Phonon read out : largest number of quanta (~10,000 phonons / keV)
 - In principle best energy resolution can be achieved in phonon channel
 - Luke-phonon readout will provide ionization energy and position information
- No further background contamination through circulation loop : no Xe circulation
- Optimal detector design for low background experiment
 - Possible container free design
 - No concerns about outgassing and leaks

Imagine

11



Short History of Solid Xenon (Argon)

12

1999@Japan

Successful contact of thin Solid Xenon to an ionization sensor using carbon graphite film






Figure 3: the first xenon layer grown on carbon graphite

2004@TAMU

Ionization readout from Solid Ar (not Xe). Failed to grow large crystals



1994@FNAL


Solid Argon detector wasn't successful

2008@FNAL

Solid Xenon Project initiated


2004@Syracuse

A student grew large Xenon crystals everyday without any problem for medical research



Solid Xenon R&D Phase-I

13




- Reproduce Syracuse setup
- Grow ~kg size of solid xenon
(inner chamber capacity : ~10kg)

2008.May : FCPA Review

2008.June : R&D phase-I approved

2008.July-Sep : Design solid xenon system
to fulfill FNAL **safety regulation**

2008.Oct - 2009.Jan : Mechanical engineering
of Stainless Steel chambers (external vendor)



Phase-I: Lab-F@Fermilab

14


Jan - Mar 2009: Install solid xenon chamber at Fermilab (Lab-F “clean room”)

April - June 2009: System check and prepare for the Fermilab **safety review**

- ➔ Top flange fabrication at FNAL mechanics shop
- ➔ Assemble parts, leak test and system check
- ➔ Grow sub-kg size xenon crystal: understand systematic issue
- ➔ **Grow ~kg size xenon crystal:** optical property check

May - June 2009: Design phase-II systems for large crystal grow and signal readout

- (1) Prepare full prescription for growing solid xenon crystal
- (2) Systematic study for growing solid argon crystal
- (3) **Propose Phase-II**



Phase-I: Xenon Chambers

15



FERMILAB CENTER FOR PARTICLE ASTROPHYSICS RETREAT, JONGHEE YOO (FERMILAB)

1. Scintillation / Ionization readout

- scintillation readout using standard photon sensors (PMT, APDs...)
- ionization readout by drifting electrons (grid mesh)
- use phase-I safety chamber (with high purity quartz vessels)
- xenon purification systems (commercial purifier or build one)



2. Demonstrate large solid xenon crystal growth (~100 kg)

- make a full prescription for growing large solid xenon
- crystal orientation measurement ?

3. Design 10 kg phase prototype detector

Fermilab Liquid Argon TPC test Facility

17




© ArTPC@FNAL



Xenon Purity Monitor

18


* Electronegative purity monitor



Material test chamber



* Optical purity monitor (new R&D)





- A new idea for optical purity monitor !
- Integrate Cavity Absorption Measurement (ICAM) → used for oceanography to measure the purity of the water
- R&D issue to find out proper reflector for xenon scintillation light ($\sim 175\text{nm}$)
- NIM paper in preparation for the R&D proposal

Phase-II: Big Vacuum Chamber

19

Used **vacuum chamber** at a mechanical vendor: perfect fit for our phase-II R&D
→ Save substantial amount of design/mechanical engineering time and budget



Glass windows




Synthetic quartz vessel
COUPP (@FNAL)

Phase-II: Orientations of the Xenon Crystal ?

20

- Use the phase-I chambers with minor modification
- X-ray / neutron diffraction measure
- Epitaxial growth of the crystal ?



E. Hurlimann et.al. (1992)

Scintillation light

PMT/QUPID/APD ...

Nano-wire / quantum dots (?)

Transition Edge Sensor (TES)

Ionization Readout

Grid mash

Electrodes / Si detector ...

TES

Phonon Readout

Cryogenics lab

- Use our collaboration (U.Florida)

Cold and warm electronics

- Well studied in CDMS collaboration

Ballistic phonon readout with TES

→ opens a whole new world!

For phonon readout

Attaching sensors on the solid xenon is the crucial step

Development of a Solid Xe Ionization Chamber

H.Nawa Y.Tamagawa M.Miyajima
Department of Applied Physics, Fukui University
9-1, Fukui 3-chome, Fukui, 910-8507, Japan




Figure 2: The temperature distribution near the contact surface of solid xenon and metal

(a).Perfect contact, (b).Imperfect contact




Figure 1: Schematic drawing for a solid xenon ionization chamber and a gas handling system




Figure 3: the first xenon layer grown on carbon graphite

Fermilab

Jonghee Yoo / Dan Bauer / Rich Schmitt / Mike Sarychev

University of Florida

Tarek Saab / Durdana Balakishieva

Texas A&M University

Rupak Mahapatra

MIT

Enectali Figueroa-Feliciano / James Kerman

Columbia

Prof. Elena Aprile (+ a postdoc and a graduate student)

Interest:

Prof. Mitsuhiro Miyajima (Waseda University)

Dr. Alexander Bolozdynya (CWRU)

Specific Questions by FCPA Directorate

23

(1) Why Fermilab?

- a) dark matter, axions, neutrino
- b) safety, infrastructure, engineers, technicians, pseudo-infinite support
- c) the idea is born at FNAL

(2) What are the risks

a) Technical

Xenon crystal will grow, but how big and how pure?

Phonon readout ? → crucial for $0\nu 2\beta$ decay experiment

b) Management/budget

Xenon price getting cheaper but still very expensive (\$4,000~\$3,000/kg)!

→ Any xenon at the Lab garage?

(3) What are the next steps?

Phase-I

Demonstrate ~kg solid xenon (hopefully within June 2009)

Write a full prescription of making xenon crystal

Phase-II

Move the all setup to PAB (currently at Lab-F)

Scintillation/ionization readout (Ar / Xe)

→ measure optical and electric characteristics

Demonstrate scalability of order 100 kg